

ONPRC Module 4B Saturday Academy

Cryopreservation & Endangered Species

| <u>Main Question</u> | <u>Laboratory Questions</u> |
|--|---|
| <p>What is “cryopreservation” and how can we use it to preserve fertility in endangered species?</p> | <ul style="list-style-type: none"> • What are the main causes for cell or tissue damages cryopreservation? • What are the 2 main techniques used in cryopreservation to minimize freezing-induced damages? • What are cryopreservation agents and how do they work? • How can we preserve fertility in endangered species using cryopreservation? |

What is cryopreservation?


Cryopreservation is a process where cells or whole tissues are preserved by cooling to low sub-zero temperatures (usually the temperature of liquid nitrogen, -196°C). At these low temperatures, any biological activity, including the biochemical reactions, that would lead to cell death is effectively stopped.



This work is licensed under a [Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License](https://creativecommons.org/licenses/by-nc-sa/4.0/).


Cryopreservation Vocabulary

Crystalline structure:



A solid composed of molecules or atoms arranged in a very specific and orderly way that displays perfect symmetry

Amorphous state (glass transition):



A solid in which there is no specific order of the atoms. Usually resulted from rapid cooling (before the molecule has enough time to organize into an orderly structure (crystal)).

<http://upload.wikimedia.org/wikipedia/commons/e/e9/Sodium-chloride-3D-ionic.png>

<http://en.wikipedia.org/wiki/Glass#mediaviewer/File:Silica.svg> Both figures - Public Domain

Viscous: Thickness of a certain fluid. Water is "thin", having a lower viscosity, while honey is "thick", having a higher viscosity.

Cryoprotective agents (CPAs): A substance that is used to protect biological tissue from freezing damage (damage due to ice formation). CPAs are usually very thick and have high viscosity.

Dehydration: The removal of water from a cell, tissue, or organism.

Liquid Nitrogen: Nitrogen in a liquid state at a very low temperature (-196°C). It can cause rapid freezing on contact with living tissue, which may lead to frostbite (use with caution! Always use goggles and thick gloves).

Vitrification: The transition of a substance into a glass (amorphous state). Cells in a solution of water and CPAs can be vitrified by rapid freezing in liquid nitrogen.

Devitrification: The growth of ice crystals during rewarming of a vitrified solution.

Supercooled liquid: Liquid at normal freezing temperature without ice formation.

Slow rate freeze: Freezing of biological samples using programmable steps at a very slow declined rate of temperature.

Seeding: Seeding is the process of inducing ice crystal formation outside the cell during slow rate freeze.



<http://www.biopoliticaltimes.org/img/original/egg%20freezing.jpg>

Cryopreserved Endangered Species



Cheetah


<http://nationalzoo.si.edu/SCBI/EndangeredSpecies/Cheetah/>

Cryobiology helped Zoo scientists breed black-footed ferrets, once thought to be extinct.

(*Mehgan Murphy/NZP*)

<http://nationalzoo.si.edu/publications/zoogoer/2010/1/FrozenAssets.cfm>





Cryopreservation of sperm and embryos is being used to preserve endangered species. They are stored in liquid nitrogen – “a zoo in a freezer”. Sperm can be thawed and used for mating with a female. Embryos can be thawed and transferred to a female for live birth.

Drawing: Joel Ito, Medical Illustrator ONPRC



Timu - baby lowland gorilla,
Photo Courtesy of Cincinnati Zoo



Baby Panda Bao Bao Photo: Abby Wood,
Smithsonian's National Zoo

<http://nationalzoo.si.edu/Animals/GiantPandas/PandaUpdates/>

Cryopreservation Laboratory Experiments

Vitrification Experiment

- Label 1 round bottom plastic tube with "0%."
- Label 4 round bottom glass tubes with "20%, 40%, 60%, 80%." Place all labeled tubes in a rack.
- Add 3ml water to a round bottom plastic tube labeled with "0%" glycerol.
- Make up 20, 40, 60, and 80% glycerol in water for a total volume of 3ml (First do the calculation (table below) and try to figure out how much glycerol and water you will need to make up each solution).

| | 0% | 20% | 40% | 60% | 80% |
|---------------------|-----|-----|-----|-----|-----|
| Glycerol (Gly) (ml) | 0.0 | | | | |
| Water (ml) | 3.0 | | | | |
| Total Volume (ml) | 3.0 | 3.0 | 3.0 | 3.0 | 3.0 |

- Mix each solution thoroughly with a transfer pipette (pipette the solution up and down 20 times).
- Set the timer for 30 seconds. Wear safety goggles and thick gloves and hold one tube (starting with the lowest %) with forceps. Submerge the tube into liquid nitrogen for 30 seconds. Make sure that all the solution is submerged below the surface of liquid nitrogen.
- At 30 seconds, carefully take the tube out of liquid nitrogen and observe the solution. Ask yourself, "is the solution vitrified?" How are you tell? Record your results in the table below.
- Repeat the last 2 steps for all your tubes.

Results

| Glycerol | 0% | 20% | 40% | 60% | 80% |
|----------------------------|----|-----|-----|-----|-----|
| Evidence of Ice? Yes or No | | | | | |
| Vitrified? Yes or No | | | | | |

Students Notes or Questions:

Cryopreservation Laboratory Experiments

Exercise

Question 1

How can you tell that a solution is vitrified successfully?

Question 2

What is the percentage of glycerol (CPA) needed for successful vitrification?

Question 3

Give an example of cryopreservation in nature.

Question 4

How can cryopreservation be used to help endangered species?

Question 5

Why is cryopreservation important for preservation of endangered species?