



NeuroQuant Scanner Parameters

V3.1, Rev 00

February 2021

Attention NeuroQuant Users:

1. NeuroQuant software requires the use of image series acquired with specified protocols as listed below.
2. Deviation from recommended protocols can result in reduced quality of NeuroQuant output.
3. Noisy scan, motion and other imaging artifacts may negatively affect results.
4. NeuroQuant requires all scan data to be in raw DICOM format. DICOM data must not be compressed in any way (e.g., no JPEG compression).
5. The software provides a data filtering mechanism, which filters scans based on parameters prior to initiating processing.
6. Do not change scan parameters.
7. MRI software updates may alter saved protocols, verify after all updates.
8. If you plan on using parallel Imaging on a 1.5T MRI scanner we will need to validate the series prior to your using it clinically. All data will be processed through NeuroQuant under research mode to avoid the standard dicom protocol checker. The raw image data noise, contrast and the segmentation results for three separate scans will be evaluated by experienced CorTechs data review team for acceptance. Please contact the support department to validate your series.

NeuroQuant / GE T1 Settings

<i>FIELD STRENGTH</i>	<i>1.5T</i>	<i>1.5T</i>	<i>3T</i>	<i>3T</i>
<i>Mode</i>	3D	*3D Parallel Imaging	3D	3D Parallel Imaging
<i>Plane</i>	Sagittal	Sagittal	Sagittal	Sagittal
<i>Pulse Sequence</i>	SPGR	SPGR	SPGR	SPGR
<i>Gradient Mode</i>	Zoom	Zoom	Zoom	Zoom
<i>Imaging Options</i>	EDR, IrP	EDR, IrP, ACC*	EDR, IrP	EDR, IrP, ACC
<i>Frequency</i>	192	192	192	192
<i>Phase</i>	192	192	192	192
<i>Frequency Direction</i>	S/I	S/I	S/I	S/I
<i>FOV</i>	24 – 25.6 cm	24 – 25.6 cm	24 – 25.6 cm	24 – 25.6 cm
<i>Phase FOV</i>	100%	100%	100%	100%
<i>Slice Thickness</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>Spacing Between Slices</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>Locs/Slab (Number of slices)</i>	160-170	160-170	160-170	160-170
<i>NEX</i>	1	1	1	1
<i>Percent Sampling</i>	100 %	100 %	100 %	100 %
<i>TE</i>	M in Full	M in Full	M in Full	M in Full
<i>Prep Time (TI)</i>	500 ms	500 ms	600 ms	600 ms
<i>Flip Angle</i>	10°	10°	8°	8°
<i>Shim</i>	Auto	Auto	Auto	Auto
<i>Bandwidth</i>	~ 16 KHz	~ 16 KHz	~ 31 KHz	~ 31 KHz

<i>Phased Array Uniformity Enhancement (PURE)</i>	On	On	On	On
<i>ASSET* Parallel Imaging Factor</i>	None	*Lowest acceptable factor	None	1.8

NeuroQuant / GE T1 Settings

Notes

1. *ASSET (ACC_GEMS) Parallel imaging compatible with Discovery 750 3T, Signa 3T, and Signa HDxt 3T scanners. NeuroQuant 2.0 or greater required
2. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For further information see page one, item 8.
3. Ensure Surface Coil Intensity Correction (SCIC) is off
4. CUBE sequence is incompatible with NeuroQuant 3D T1 sequence
5. Do not add post processing filters
6. Do not use Zip512

NeuroQuant / Philips T1 Settings

<i>FIELD STRENGTH</i>	<i>1.5T</i>	<i>*1.5T</i>	<i>3T</i>	<i>3T</i>
<i>Mode</i>	3D	3D Parallel Imaging	3D	3D Parallel Imaging
<i>Plane</i>	Sagittal	Sagittal	Sagittal	Sagittal
<i>Technique</i>	FFE	FFE	FFE	FFE
<i>Contrast Enhancement</i>	T1	T1	T1	T1
<i>Fast Imaging Mode</i>	TFE	TFE	TFE	TFE
<i>FOV</i>	240-256 mm	240-256 mm	240-256 mm	240-256 mm
<i>RFOV</i>	100%	100%	100%	100%
<i>Matrix Scan</i>	192	192	192	192
<i>Stacks</i>	1	1	1	1
<i>Scan Percentage</i>	100%	100%	100%	100%
<i>Slices</i>	160-170	160-170	160-170	160-170
<i>Slice Thickness</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>Spacing Between Slices</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>TFE Shot Interval</i>	2300 ms	2300 ms	2500 ms	2500 ms
<i>Echoes</i>	1	1	1	1
<i>TE</i>	4 ms	4 ms	Shortest	Shortest
<i>TR</i>	Shortest	Shortest	Shortest	Shortest
<i>Flip Angle</i>	8°	8°	9°	9°
<i>TFE Prepulse Delay / TI</i>	1000 ms	1000 ms	900 ms	900 ms
<i>Water Fat shift</i>	Maximum	Maximum	Maximum	Maximum
<i>Gradient Mode</i>	Default	Default	Maximum	Maximum
<i>*SENSE Parallel Acceleration Factor</i>	None	*Lowest acceptable factor	None	1.8

NeuroQuant / Philips T1 Settings

Notes

1. 1.5T scan protocol is based on Philips Achieva scanner with software version 3.2.2
2. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For further information see page one, item 8.
3. SENSE Parallel imaging compatible with Achieva, Ingenia, and Intera **3T** scanners
4. CLEAR (Homogeneity Correction) must be **ON**
5. DICOM connection must be Classic
6. Do not add post processing filters

NeuroQuant / Siemens T1 Settings

<i>FIELD STRENGTH</i>	<i>1.5T</i>	<i>1.5T</i>	<i>3T</i>	<i>3T</i>
<i>Mode</i>	3D	3D Parallel Imaging	3D	3D Parallel Imaging
<i>Plane</i>	Sagittal	Sagittal	Sagittal	Sagittal
<i>Sequence</i>	MPRAGE	MPRAGE	MPRAGE	MPRAGE
<i>FOV Phase</i>	100%	100%	100%	100%
<i>FOV Read</i>	240-256 mm	240-256 mm	240-256 mm	240-256 mm
<i>Base Resolution</i>	192	192	192	192
<i>Slice Thickness</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>Spacing Between Slices</i>	1.2 mm	1.2 mm	1.2 mm	1.2 mm
<i>Phase Oversampling</i>	0%	0%	0%	0%
<i>Slice Oversampling</i>	0%	0%	0%	0%
<i>Slices per Slab</i>	160-170	160-170	160-170	160-170
<i>Averages</i>	1	1	1	1
<i>Concentrations</i>	1	1	1	1
<i>Gradient Mode</i>	Fast	Fast	Fast	Fast
<i>Mag. Preparation</i>	Non-Sel. IR	Non-Sel. IR	Non-Sel. IR	Non-Sel. IR
<i>RF Pulse Type</i>	Fast	Fast	Fast	Fast
<i>Excitation</i>	Non-Sel. IR	Non-Sel. IR	Non-Sel. IR	Non-Sel. IR
<i>TR</i>	2400 ms	2400 ms	2300 ms	2300 ms
<i>TE</i>	Minimum	Minimum	Minimum	Minimum
<i>TI</i>	1000 ms	1000 ms	900 ms	900 ms
<i>Flip Angle</i>	8°	8°	9°	9°
<i>Bandwidth</i>	180 HZ/Pix	180 HZ/Pix	240 HZ/Pix	240 HZ/Pix
<i>Filter</i>	Non	Non	Non	Non

<i>Parallel Technique iPAT</i>	Off	Validation only	Off	On
<i>Max Parallel Acceleration Factor</i>	N/A	*Lowest acceptable factor	N/A	GRAPPA 2

NeuroQuant / Siemens T1 Settings

Notes

1. GRAPPA parallel imaging compatible with Verio, Trio Tim and Skyra 3T scanners only
2. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For further information see page one, item 8.
3. Pre-scan Normalize must be **ON**
4. Do not add post processing filters

NeuroQuant / Canon (Toshiba) T1 Settings

<i>FIELD STRENGTH</i>	<i>*1.5T</i>	<i>**3T</i>
<i>Mode</i>	3D	3D
<i>Plane</i>	Sagittal	Sagittal
<i>Sequence Option</i>	FFE3D	FFE3D
<i>Sequence Option</i>	Basic	Default
<i>Matrix</i>	192x192	192x192
<i>FOV</i>	24 cm	24 cm
<i>Slice Thickness</i>	1.2 mm	1.2 mm
<i>Spacing Between Slices</i>	1.2 mm	1.2 mm
<i>RF Type</i>	Fine	Normal
<i>Gradient Type</i>	Normal	Normal
<i>TE</i>	2.9 ms	3.2 ms
<i>TR</i>	7.1 ms	7.3 ms
<i>IR Pulse</i>	Non-Sel.	Non-Sel.
<i>TI</i>	800 ms	900 ms
<i>Flip Angle</i>	9°	9°
<i>Bandwidth</i>	217 HZ/pix	244 HZ/pix
<i>No Wrap</i>	1	1
<i>Speeder</i>	PE 1	PE 1
<i>Number of Segments</i>	2 / Sequential	1
<i>Type of Shot Interval</i>	1600 ms	2300 ms
<i>IDC (Gradient Correction)</i>	VMC	VMC
<i>NAQ</i>	1	1

NeuroQuant / Canon (Toshiba) T1 Settings

Notes

1. Settings are for Toshiba Titan MRI scanners
2. Do not add post processing filters

*1.5T settings only to be used with MPower v3.6 software or later. NeuroQuant 2.1 or greater required.

**3T settings only to be used with Mpower v3.5 software or later. NeuroQuant 2.0 or greater required.

NeuroQuant / Hitachi T1 Settings

<i>FIELD STRENGTH</i>	1.2T	1.5T	3T
<i>Acquisition Type</i>	3D	3D	3D
<i>Plane</i>	Sagittal	Sagittal	Sagittal
<i>Mode</i>	GEIR	GEIR	GEIR
<i>Sequence</i>	RSSG	RSSG	RSSG
<i>FOV</i>	240 - 256 mm	240 - 256 mm	240 - 256 mm
<i>Acquisition Matrix</i>	192 x 192	192 x 192	192 x 192
<i>Slice Thickness</i>	1.2 mm	1.2 mm	1.2 mm
<i>Spacing Between Slices</i>	1.2 mm	1.2 mm	1.2 mm
<i>NSA</i>	1	1	1
<i>Slices Per Slab (# of Slices)</i>	160 - 170	160 - 170	160 - 170
<i>TR</i>	8.7 ms	8.6 ms	7.2 ms
<i>TE</i>	4 ms	4 ms	2.9 ms
<i>TI</i>	550 ms	550 ms	900 ms
<i>Flip Angle</i>	12°	12°	10°
<i>Bandwidth</i>	53.2 KHz	58.2 KHz	92.6 KHz
<i>Parallel Imaging RAPID Phase & Slice</i>	1, 1	1, 1	1, 1
<i>Filter Option VIVID</i>	3	1	1

Notes

1. Gradient Correction should be **ON** for all magnetic field strengths & modes
2. Do not use post processing filters other than Vivid



NeuroQuant MS Scanner Parameters

V3.1, Rev 00

February 2021

Attention NeuroQuant MS Users:

Prior to using NeuroQuant MS module, the same version of NeuroQuant (3D T1) and NeuroQuant MS (T2 FLAIR) scanner settings must be configured on your scanner(s).

NOTE: You do not need both the 2D and the 3D T2 FLAIR sequences to generate NeuroQuant MS reports. You need the NeuroQuant Sagittal 3D T1 and either the 2D or the 3D T2 FLAIR sequences.

If you plan on using parallel Imaging on a 1.5T MRI scanner we will need to validate the series prior to your using it clinically. All data will be processed through NeuroQuant under research mode to avoid the standard dicom protocol checker. The raw image data noise, contrast and the segmentation results for three separate scans will be evaluated by experienced CorTechs data review team for acceptance. Please contact the support department to validate your series.

NeuroQuant MS / GE T2 FLAIR Settings

<i>FIELD STRENGTH</i>	1.5T		3T	
<i>Mode</i>	2D	3D	2D	3D
<i>Plane</i>	Axial	Sagittal	Axial	Sagittal
<i>Frequency</i>	256	256	256	256
<i>Phase</i>	256	256	256	256
<i>FOV</i>	220 - 240	220 - 240	220 - 240	220 - 240
<i>Slice Thickness</i>	≤ 3 mm	1.2 mm	< 3 mm	1.2 mm
<i>TE</i>	127 ms	90 ms	147 ms	130 ms
<i>TR</i>	8,000 ms	6,900 ms	11,000 ms	6,000 ms
<i>Flip Angle</i>	90°	90°	90°	90°
<i>TI</i>	2,000 ms	2,030 ms	2,250 ms	1,854
<i>Bandwidth</i>	195 Hz/px or 50 KHz	122 Hz/px or 31 KHz	122 Hz/px or 31 KHz	122 Hz/px or 31 KHz
<i>ASSET Parallel Imaging Factor</i>	NONE	*Lowest factor acceptable	2	2

Notes

1. ASSET parallel imaging compatible with Discovery 750 3T, Signa 3T, and Signa HDxt 3T scanners only
2. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For further information see page one, item 8.
3. Variable flip angle flag should be off
4. Variable Bandwidth should be off
5. Do not use post processing filters
6. ZIP512 not supported

NeuroQuant MS / Philips T2 FLAIR Settings

<i>FIELD STRENGTH</i>	1.5T		3T	
<i>Mode</i>	2D	3D	2D	3D
<i>Plane</i>	Axial	Sagittal	Axial	Sagittal
<i>Gradient Mode</i>	Default	Default	Default	Default
<i>Fast Imaging Mode</i>	Default	Default	FSE	FSE
<i>In Plane Pixel Resolution</i>	1.0 mm x 1.0 mm	1.0 mm x 1.0 mm	1.0 mm x 1.0 mm	1.0 mm x 1.0 mm
<i>Percent Phase FOV</i>	100%	100%	100%	100%
<i>FOV</i>	22 – 25.6 cm	22 – 25.6 cm	22 – 25.6 cm	22 – 25.6 cm
<i>Slice Thickness</i>	≤ 3 mm	2.4 mm	≤ 3 mm	1.2 mm
<i>TE</i>	120 ms	296 ms	90 ms	274 ms
<i>TR</i>	6,000 ms	8,000 ms	9,000 ms	4,800 ms
<i>Inversion Time/TI</i>	2,000 ms	2,400 ms	2,500 ms	1,650 ms
<i>Flip Angle</i>	90°	90°	90°	90°
<i>SENSE Parallel Imaging Factor</i>	1	*Lowest factor acceptable	1.6	≤ 2

Notes

1. 1.5T scan protocol is based on Philips Achieva scanner with software version 3.2.2
2. SENSE Parallel imaging compatible with Achieva, Ingenia 3T, and Intera 3T scanners only
3. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For further information see page one, item
4. Do not use post processing filters
5. ZIP512 not supported

NeuroQuant MS / Siemens T2 FLAIR Settings

<i>FIELD STRENGTH</i>	1.5T		3T	
<i>Mode</i>	2D	3D	2D	3D
<i>Plane</i>	Axial	Sagittal	Axial	Sagittal
<i>Matrix Scan</i>	256	256	256	256
<i>Base Resolution</i>	256	256	256	256
<i>Phase Resolution</i>	100%	100%	100%	100%
<i>FOV</i>	22 - 25.6	22 - 25.6	22 - 25.6	22 - 25.6
<i>Slice Thickness</i>	≤ 3 mm	1.2 mm	< 3 mm	1.2 mm
<i>Averages</i>	1	1	1	1
<i>Magn. Preparation</i>	Slice-sel IR	Slice-sel IR	Slice-sel IR	Slice-sel IR
<i>TE</i>	111 ms	358 ms	90 ms	392 ms
<i>TR</i>	9,000 ms	5,000 ms	9,000 ms	5,000 ms
<i>Inversion Time/TI</i>	2,500 ms	1,800 ms	2,500 ms	1,800 ms
<i>Flip Angle</i>	150°	120°	150°	120°
<i>Bandwidth</i>	130 Hz/px or 33 KHz	592 Hz/px or 152 KHz	222 Hz/px or 57 KHz	780 Hz/px or 32 KHz
<i>GRAPPA Parallel Imaging Factor</i>	NONE	*Lowest acceptable factor	2	2

Notes

1. GRAPPA parallel imaging compatible with Verio 3T, Trio Tim 3T and Skyra 3T scanners only
2. *In order to use parallel imaging with a 1.5T MRI sequence, it must be validated by the Cortechs team. Lower PI factors produce better results. For more information see page one, item 8.
3. Do not use post processing filters
4. ZIP512 not supported

NeuroQuant MS / Canon (Toshiba) T2 FLAIR Settings

FIELD STRENGTH	1.5T		3T	
Mode	2D	3D	2D	3D
Plane	Axial	Sagittal	Axial	Sagittal
Sequence	FSE2D+13	FASE3D+5.5	FSE2D+12	FASE3D+5.5
Sequence Option	Basic	MPV T2	Basic	MPV T2
In Plane Pixel Resolution	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm
FOV	22 - 25.6 cm	22 - 25.6 cm	22 - 25.6 cm	22 - 25.6 cm
Slice Thickness	3 mm	1.2 mm	3 mm	1.2 mm
Gradient Type	SLT	N/A	Normal	N/A
RF Type	Normal	N/A	Low SAR	N/A
TE	104 ms	352 ms	96 ms	401.5 ms
TR	10,000 ms	7,000 ms	9,000 ms	6,000 ms
TI	2,500 ms	2,200 ms	2,500ms	2,000 ms
Flip Angle	90°	90°	90°	90°
FLOP	180°	180°	140°	180°
IR Pulse	Robust	Robust	Standard	Robust
No Wrap	1	SE 1.2	None	SE 1.2
Speeder	PE 1	PE 2 SE 2	PE 2	PE 2 SE 2
Number of Shots	N/A	1	N/A	1
Bandwidth	139.5 Hz/pix	651 Hz/pix	195.3 Hz/pix	651 Hz/pix

Notes

1. Fat Suppression should be off
2. IDC (Gradient Correction) – set to VMC
3. Do not use post processing filters
4. ZIP512 not supported

NeuroQuant MS / Hitachi T2 FLAIR Settings

<i>FIELD STRENGTH</i>	1.2T		1.5T		3T	
<i>Mode</i>	2D	3D	2D	3D	2D	3D
<i>Plane</i>	Axial	Sagittal	Axial	Sagittal	Axial	Sagittal
<i>Sequence</i>	opFIR	isoFSE	opFIR	isoFSE	FIR	isoFSE
<i>In Plane Pixel Resolution</i>	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm	1 mm x 1 mm
<i>FOV</i>	22 cm	24 - 26 cm	22 cm	24- 26 cm	22 cm	24 - 26 cm
<i>Slice Thickness</i>	≤ 3 mm	1.5 acquired / 1 recon	≤ 3 mm	1.2 acquired / 1 recon	≤ 3 mm	1.2 acquired / 1 recon
<i>TE</i>	132 ms	390 ms	132 ms	390 ms	116 ms	438 ms
<i>TR</i>	9,130 ms	5,500 ms	9,660 ms	6,000 ms	10,000 ms	10,000 ms
<i>TI</i>	2,400 ms	1,800 ms	2,500 ms	1,860 ms	2,300 ms	2,600 ms
<i>Flip Angle</i>	90°	90°	90°	90°	90°	90°
<i>Flip Rapid Phase</i>	1	2	1.2	2	1.1	2
<i>Rapid Slice</i>	N/A	1.4	N/A	1.7	N/A	2.5
<i>Rapid Mode</i>	N/A	RCM	RCM	RCM	RCM	RCM
<i>Vivid Filter</i>	3	1	3	2	3	1
<i>Bandwidth</i>	56 KHz	96.2 KHz	65 KHz	91 KHz	80 KHz	108.7 KHz

Notes

1. Percent Phase FOV should be 100%
2. Gradient Correction should be On
3. Do not use post processing filters other than Vivid
4. ZIP512 not supported