

# Bruker Spatial Biology: NanoString nCounter Technology

@ Oregon Health & Science University

---

**Bailey Longoni**  
nCounter Regional Account Manager  
October 15<sup>th</sup> 2024



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

# We're here to help!

- **NanoString Support Group**

- Bailey Longoni – nCounter Regional Account Manager
  - [Bailey.Longoni@bruker.com](mailto:Bailey.Longoni@bruker.com)
- Wes Heydeck – Field Application Scientist
  - [Westley.Heydeck@bruker.com](mailto:Westley.Heydeck@bruker.com)
- Free online videos and training: NanoString University (<https://university.nanostring.com>)
- Support: [support.spatial@bruker.com](mailto:support.spatial@bruker.com)

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Tools for Biomarker Discovery and Translational Research

Founded: 2003  
Headquartered: Seattle, WA



> 5,000 Publications To Date



**nCounter®  
Analysis System**  
Launched 2008

**GeoMx® Digital  
Spatial Profiler**  
Launched 2019



**CosMx® Spatial  
Molecular Imager**  
Launched 2022


FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

# NanoString nCounter




# Gene expression measurement is fundamental



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

# Challenges of Current Gene Expression Workflows



qPCR, microarray & RNA-Seq workflows contain numerous steps where variability can be introduced into the data

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

A BRUKER COMPANY

# nCounter Helps RNA-Seq Users Validate Their Findings During Translational Research


Hypothesis formation using whole transcriptome discovery



Validate findings using a technology well suited for translational research

If you can't orthogonally validate your results, are they valid?

## Complex Workflow



## Sources of Variation

Technique Bias

Manual Pipetting Errors

Reverse Transcriptase Efficiency

Amplification Bias/Errors

## RNAseq Struggles with Clinical Samples

### RESEARCH ARTICLE

A Comparison of RNA-Seq Results from Paired Formalin-Fixed Paraffin-Embedded and Fresh-Frozen Glioblastoma Tissue Samples

Anna Esteve-Codina<sup>1</sup>, Oriol Arpi<sup>2</sup>, Maria Martinez-García<sup>3</sup>, Estela Pineda<sup>4</sup>, Mar Mallo<sup>5</sup>, Marta Gut<sup>6</sup>, Cristina Carrato<sup>7</sup>, Anna Rovira<sup>2</sup>, Raquel Lopez<sup>8</sup>, Avelina Tortosa<sup>9</sup>, Marc Dabad<sup>1</sup>, Sonia Del Barco<sup>10</sup>, Simon Heath<sup>1</sup>, Silvia Bague<sup>11</sup>, Teresa Ribalta<sup>12</sup>, Francesc Alameda<sup>13</sup>, Nuria de la Iglesia<sup>14</sup>, Carmen Balaña<sup>15\*</sup>, on behalf of the GLIOCAT Group<sup>11</sup>


In this RNA-Seq study, informative results were obtained from only 3 out of 11 matched FF& FFPE samples

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

# Molecules That Count<sup>®</sup>

Detection of individual nucleic acid molecules via hybridization to a fluorescent barcode







FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY




# Standard nCounter barcode chemistry



Barcode	Identity
	XLSA
	FOX5
	PDCD1


FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Standard nCounter barcode chemistry



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# The nCounter<sup>®</sup> miRNA Assay Adapts the Dual Probe System to Detect Small RNAs




FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

# Automated digital profiling in two simple steps

 **5** min  
HANDS-ON

Day **1**

 **10** min  
HANDS-ON

Day **2**  
AUTOMATED

## Hybridize



- Flexible sample requirements
- Only 4 pipetting steps
- No amplification
- 800 hybridizations in single tube

## Purify & Count

*nCounter<sup>®</sup> SPRINT*



- Sensitive
- Precise
- Quantitative
- Simple

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.


**nanoString**  
A BRUKER COMPANY

nanosString<sup>TM</sup>

The logo features the text "nanosString" in a black, sans-serif font. A horizontal line is positioned below the text, starting from the left and ending at the right edge of the "g". Five small, light green spheres are placed on this line, evenly spaced between the "n", "a", "n", "o", and "S". A small "TM" trademark symbol is located at the bottom right of the "g".



# Digital output: 1 molecule = 1 count



Codes are Counted and Tabulated

Barcode	Counts	Identity
	3	XLSA
	2	FOX5
	1	INSULIN



	Sample 1	Sample 2	Sample 3
Assay number	1	2	3
Sample ID	PW1	PW2	PW3
Reporter Counts			
Gene Name			
POS_A(128)	16874	16930	17013
POS_B(32)	4122	4134	4105
POS_C(8)	1144	1087	1078
POS_D(2)	379	362	357
POS_E(0.5)	82	85	55
POS_F(0.125)	36	38	28
NEG_A(0)	14	5	5
NEG_B(0)	15	17	16
NEG_C(0)	14	13	7
NEG_D(0)	23	20	31
NEG_E(0)	20	13	18
NEG_F(0)	13	7	8
NEG_G(0)	16	9	11
NEG_H(0)	27	20	9
G6PD	9958	11129	11059
GUSB	2443	2747	2843
HPRT1	3801	4311	4184
TBP	1221	1311	1359
POLR1B	1068	1124	1114

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Data Analysis


# Data Analysis Analysis Tools Eliminate Time Spent Waiting for Results

## nSolver



Differential expression analysis, PCA analysis, cell type profiling, pathway scoring, etc.

## ROSALIND®: Cloud-Based Unified Data Analysis Platform



Evolving knowledge base (50+)  
Shared collaboration spaces




# Automated Data Analysis with nSolver™ Software


- Rapid data QC and normalization
- Advanced analysis modules look at differential expression, pathways, cell type profiling and more
- Compatible with standard analysis programs for existing workflow





# Digital Data Output: Simple Visualizations via nSolver™ 4.0




Box Plot by Treatment Group



Scatter Plot of All Genes



New Universal Advanced Analysis Module



Clustering and Heat Map


FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.




# Simplified Data Analysis – nSolver™ Advanced Analysis

nCounter Advanced Analysis is a free, **easy-to-use** add-on to nSolver Basic Software for deeper data insights based on robust R statistics.


## Differential Expression Analysis




## PCA Analysis



## Cell Type Profiling



## Pathway Scoring




FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.



# INTERACTIVE VISUALIZATIONS

## INTERPRET YOUR DATA LIKE A BIOINFORMATICIAN



Empower researchers to:

- ✓ Reduce the time required to get from sequencing to data interpretation
- ✓ Enjoy an engaging user experience to interactively explore results and pathways
- ✓ Setup experiments and explore interactive results the very same day
- ✓ Create new comparisons and change cut-offs in minutes
- ✓ No programming expertise required



# COLLABORATE EFFECTIVELY

## AGGREGATE AND ANALYZE SILOED DATASETS



- Breakdown Barriers to Collaboration
  - ✓ Collaborate in real-time, any where in the world, with consistent access to shared experiments for all participants
  - ✓ No need to download data or transfer data
  - ✓ Collaborate with an engaging user experience, not offline CSV files

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.



**nanoString**  
A BRUKER COMPANY

# Performance

# Technical reproducibility





qPCR



\*CNRQ=Calibrated Normalized Relative Quantity (qbase+ software)  
Data shown in this experiment are representative of many runs

nCounter



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.









# Reproducibility between researchers



qPCR



nCounter




FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.



# Reproducibility between researchers





qPCR



**No need for technical replicates!**

nCounter




FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.



# High sensitivity and dynamic range




## Sensitivity, Dynamic Range, Reproducibility



- Reproducible and linear over a wide dynamic range (5-6 logs)
- Note the reproducibility even at very low expression levels

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Consistent performance with varying sample preparation




FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Sample Inputs




# Challenging Sample Types Are No Challenge for nCounter

## Broad Sample Compatibility



## Even Low-Quality Samples Succeed



*Decades Old FFPE Directly Correlates with FF*

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# General guidelines for NanoString nCounter samples


- Good quality RNA
  - $A_{260}/A_{280}$  ratios ~ 2.0
  - $A_{260}/A_{230}$  ratios ~ 2.0
  - DV200 = 50%
- Typical inputs for gene expression assays from high quality samples
  - SPRINT system: 50 ng (25 - 150 ng)
- Keep in mind only ~5% of total RNA is mRNA!



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY


# RNA is RNA is RNA...It's not just human



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Gene Expression Panels

# Panels – Ready to use – or make your own!




# And more!



**nanoString**  
A BRUKER COMPANY



# Cardiovascular Disease Pathophysiology Panel




Cardiovascular Pathology	Cardiovascular Physiology	Vascular Inflammation	Cellular Aging & Renewal	Metabolism	Mechano Signaling	Regulatory Signaling	Epigenetic Remodeling
Atherosclerosis Cardiac Hypertrophy Cardiomyopathy Foam Cell Formation Ischemia Myocarditis Thrombosis Hemostasis Endocarditis Pericarditis Hypertension	Angiotensin System Cardiac Muscle Contraction Cardiac Electrophysiology GABAergic Signaling Vascular Smooth Muscle Contraction Vasopressin System Cardiac Morphogenesis ER Stress	eNOS Activation IL-1 Signaling IL-6 Signaling Other Cytokine Signaling Immune Cell Infiltration JAK-STAT Signaling mTOR Signaling NF-kappaB Signaling PI3K-AKT Signaling PPAR Signaling TLR Signaling TNF Signaling Checkpoint Signaling	Apoptosis Autophagy Cell Cycle Senescence & Quiescence Telomere Maintenance	Fatty Acid Metabolism Glucose Metabolism Hypoxia Response Lipid Metabolism Cholesterol Metabolism Lipoprotein Clearance Oxidative Stress Response	ECM Remodeling Hippo Signaling Integrin Signaling Rho ROCK Signaling	Calcium Signaling EGFR Signaling MAPK Signaling Notch Signaling TGF-beta Signaling VEGF Signaling Wnt Signaling	Histone Modifications Acetyl Transferases Deacetylases Methyl Transferases

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.



# Cardiovascular Disease Pathophysiology Panel



Cardiovascular Pathology	Cardiovascular Physiology	Vascular Inflammation	Cellular Aging & Renewal	Metabolism	Mechano Signaling	Regulatory Signaling	Epigenetic Remodeling
Atherosclerosis Cardiac Hypertrophy Cardiomyopathy Foam Cell Formation Ischemia Myocarditis Thrombosis Hemostasis Endocarditis Pericarditis Hypertension	Angiotensin System Contraction Cardiac Electrophysiology GABAergic Signaling Vascular Smooth Muscle Contraction Vasopressin System Cardiac Morphogenesis ER Stress	NOS Activation IL-1 Signaling IL-6 Signaling Other Cytokine Signaling Immune Cell Infiltration JAK-STAT Signaling mTOR Signaling NF-kappaB Signaling PI3K-AKT Signaling PPAR Signaling TLR Signaling TNF Signaling Checkpoint Signaling	Apoptosis Autophagy Cell Cycle Senescence & Quiescence Telomere Maintenance	Fatty Acid Metabolism Glucose Metabolism Hypoxia Response Lipid Metabolism Cholesterol Metabolism Lipoprotein Clearance Oxidative Stress Response	ECM Remodeling Hippo Signaling Integrin Signaling Rho ROCK Signaling	Calcium Signaling EGFR Signaling MAPK Signaling Notch Signaling TGF-beta Signaling VEGF Signaling Wnt Signaling	Histone Modifications Acetyl Transferases Deacetylases Methyl Transferases

Gene	Cell Type	Internal Reference Gene	Angiotensin System	Apoptosis	Atherosclerosis
CSK		-	-	-	-
CSNK2B		-	-	-	-
CST3		-	-	-	-
CTF1		-	-	-	-
CTLA4		-	-	-	-
CTNNB1		-	-	-	-
CTNND1		-	-	-	-
CTSA		-	+	-	-
CTSD		-	+	+	-
CTSG		-	+	-	-
CTSW	Cytotoxic Cell	-	-	+	-
CTSZ		-	+	+	-
CX3CL1		-	-	-	-
CXCL12		-	-	-	-
CXCL8		-	-	-	+
CXXC1		-	-	-	-
CYBA		-	-	-	+
CYC1		-	-	-	-
CYCS		-	-	+	+
CYGB		-	-	-	-
CYP27A1		-	-	-	-
CYP2U1		-	-	-	-
CYSLTR2		-	-	-	-
DCN	Fibroblasts	-	-	-	-
DCXR		-	-	-	-
DDAH1		-	-	-	-
DDAH2		-	-	-	-

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# nCounter Panel Selection Tool

PANEL PRO TUTORIAL

## Browse Panels

[View all panels](#)

Enter panel name, gene symbol, annotation

### By application

Neuroscience  
Oncology  
Immunology

[View all >](#)

### By pathway

A1 Astrocyte  
A2 Astrocyte  
Activated

[View all >](#)

### By biological process

actin cytoskeleton organization  
actin cytoskeleton reorganization  
actin filament bundle assembly

[View all >](#)

## Search By Gene List

[Download template](#)

[Upload list](#)

Enter gene symbols, separated by commas

Values entered: 0

## Compare Panels

Enter up to 6 panels for comparison

# Panel Selection Tool

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Customize an existing panel or create a new one

## Panel Plus:

Customizable Add-on to a Panel


OR

Build Content that is Customized  
for Any Gene, Any Organism

PanCancer IO 360™  
Gene Expression Panel



Up to 55  
User-defined Targets



Custom  
CodeSets

## Turnkey Solution for Any Project

- Maximum target number: 800
- Study novel and non-model organisms
- Address transcript variants and transgenes

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

**nanoString**  
A BRUKER COMPANY

## Welcome to NanoString University



### Learn More About NanoString's Platforms

Click here to get started with introductory level courses and hear from colleagues in the field!



### Training Videos

Access all our how-to videos and full length training courses here!



### Document Library

Find nCounter and GeoMx user manuals, protocols, guidelines, and more.



### GeoMx<sup>®</sup> DSP Site Readiness Portal

Congratulations on purchasing your new GeoMx<sup>®</sup> DSP! Visit this page to start preparing for your upcoming training.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Support for Data Analysis on NanoString University



## nSolver 101: Basic Analysis

Sign up for this course to learn the basics of nSolver data analysis.


1 hr 8 min



## nSolver 201: Advanced Analysis

Learn how to take your nSolver analysis to the next level!

48 min



nanoString

Registered

Analyzing  
nCounter®  
Data with  
ROSALIND®

## Using ROSALIND® to Analyze nCounter® Data

Sign up for this course to learn how to analyze your data using the ROSALIND® cloud-based software!

1 hr 34 min

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

# Thank You!

nanoString



**15% discount on all off-the-shelf nCounter Panel orders through December 13<sup>th</sup>, 2024**

Schedule a project consultation through [bailey.longoni@bruker.com](mailto:bailey.longoni@bruker.com)

# Classifying AML Tumor Cell State to Identify Responsiveness to Venetoclax-Based Therapy

Steve Kurtz

Chris Eide

Alexandria Miller

Sud Anand

Dan Bottomly

Jeff Tyner




KNIGHT  
CANCER  
*Institute*

10/15/2024



# Acute Myeloid Leukemia (AML) results from the expansion of immature white blood cells in the bone marrow and peripheral blood



Hematopoietic stem cells are developmentally blocked and proliferate as immature myeloid cells (blasts)

Genetically heterogeneous: numerous mutations and translocations

Conventional treatment: Intensive chemotherapy of cytarabine+daunorubicin (7+3)


More recently: Hypomethylating agents (azacitidine, decitabine), Bcl2 inhibitor (venetoclax)

Survival remains poor as patients relapse

# Advances and limitations of venetoclax-based therapy for AML

Azacytidine plus the BCL2 inhibitor venetoclax (Ven) is standard-of-care for elderly or unfit patients with AML


*DiNardo et al., NEJM 2020*



Ven+Azacytidine improves response rates with certain genetic mutations (e.g. *IDH1/2*, *NPM1*), with poorer responses in *TP53* and *RAS* mutated AML

*Stahl et al., Blood Adv 2021*

*DiNardo et al., Blood 2020*



Differential sensitivity to Ven +/- Azacytidine associates with primitive but not differentiated AML tumor cell state

*Pei et al., Cancer Discov 2020*

*Zhang et al., Nat Cancer 2020*

*Kuusanmäki et al., Haematologica 2020*

# Classifying AML by tumor cell state

Single-cell RNAseq and genotyping methodology identified gene sets that distinguish cell hierarchy lineages and differentiation




*van Galen et al., Cell 2019*

Application of these signatures to AML patient samples determined tumor cell state (e.g. primitive or monocytic) associations with Ven ex vivo sensitivity

*Bottomly et al., Cancer Cell 2022*

*Zeng et al., Nature Medicine 2022*

*Eide and Kurtz et al., Blood Cancer Discov 2023*



RNAseq data

# Translating the findings: Assays to rapidly profile tumor cell state

~50-60% of newly diagnosed AML patients are elderly or unfit for intensive chemotherapy;  
these patients receive Ven-based treatment

*Latagliata et al., Annals of Onc 2006*

*Oyogoa et al., Cancers (Basel) 2023*

~30% of elderly or unfit AML patients treated with frontline Aza+Ven do not achieve complete response



*DiNardo et al., NEJM 2020*

RNAseq is gold standard for transcriptomic profiling, but analysis time/cost effectiveness of prospective single sample runs not compatible with clinical decision-making timelines in AML

Goal: Develop a rapid assay to profile cell state gene expression signatures in AML patients and align with prediction of Aza+Ven sensitivity

# Technology for Rapid Determination of Cell State

## Nanostring Technology



Designed custom Nanostring panel targeting:


AML hematopoietic cell state signatures (6 sets of 30 genes/set from van Galen et al., *Cell* 2019)

AML LSC17 “stemness” signature (17-gene set from Ng et al., *Nature* 2016)

Control genes (n=7)

# Classification of AML Tumors by Cell State Profiling

Pilot study: RNAs from 5 samples previously defined as Progenitor-like and 5 samples previously defined as Monocytic-like by deconvolution from bulk RNAseq



Samples are clustered in rows and genes are clustered, grouped by gene set. Shading reflects the z-scores for each sample for a given gene with increased (orange) or decreased (blue) color normalized gene expression.

# Design Considerations and Performance

Sample Expression Distribution is consistent across independent batch runs

QC criteria were determined based on Nanostring recommendations, existing CLIA lab protocols, and collaborator input.

Field of view (FOV) > 75%


Binding density > 0.5

Positive control linearity R2 values at least 0.95

Lower limit of detection (LOD) threshold value is 2 SD above the geometric mean of the negative controls


The lowest abundance positive control probe (0.5 fM) > LOD threshold value

## Distribution of Raw Counts Across Batches





# Nanostring Profiling on 49 RNAs with matching RNAseq data




# Nanostring Assay concordance with RNAseq

Gene set enrichment was calculated using a modified single-sample method (Singscore)

Foroutan et al., *BMC Bioinformatics* 2018

Nanostring-based results correlate with cell state scores determined from deconvolution of bulk RNAseq data

Bottomly et al., *Cancer Cell* 2022




# Changes in Tumor cell state profile on Intensive Chemotherapy: Patient 2314


At diagnosis, patient with MDS, treated with azacytidine, but transformed to AML with monocytic differentiation.

Treated with 7+3 induction, did not achieve a complete response.


Tumor profile shifted to primitive cell state following treatment.



Classifying cell state at AML diagnosis provides context for anticipating response to Aza+Ven




**Clinical AZA+VEN outcome:**  
Complete Response




**Clinical AZA+VEN outcome:**  
Refractory

## Changes in Tumor cell state profile on AZA+VEN : Patient 6282



Tumor profile shifted from primitive to monocytic cell state following treatment.

# Translational Goal





*University Shared Resources*

# NanoString Services

RNA & DNA Services of the Integrated Genomics Laboratory (IGL)



Jinah Kim, PhD  
Core Scientist  
Oct, 15<sup>th</sup>, 2024




# Advantages of NanoString Technology

- ❖ Direct quantification of RNA and miRNA using molecular barcode without the need for amplification
- ❖ Each multiplex panel can profile up to 770 genes, with the option to customize by adding up to 55 additional unique targets
- ❖ Faster than NGS (in terms of setup and processing), more sensitive than targeted RNAseq (especially for genes with low expression levels), and simpler than qPCR (as it requires no amplification and uses a direct hybridization method).

We are located in Richard T. Jones Hall (RJH), 5th floor, Room 5330.

RJH



Integrated Genomics Laboratory

# NanoString nCounter SPRINT Workstation (RJH, 5330)



**Day1:** CodeSet Hybridization using the thermocycler (**16 hours**)

**Day 2:** Sample loading into the Cartridge and running the nCounter SPRINT Profiler (**6 hours**)



# Training is required

- **Training includes:**
  - Technical Information and Background Reading
  - Virtual Meeting (45 minutes)
  - Hands-on, in-person training (45 minutes)
- **Current Training Fee:**
  - \$150 (Internal)
  - \$180 (External Academic)
  - Max 2 people per lab per training session
  - Slated to increase in FY25

**Post Training Support** - For your first NanoString run, IGL staff will shadow trainees to provide support with CodeSet hybridization, cartridge loading, and the instrument operation. Staff remain available during subsequent runs to answer questions.

# Reservation


- **Bookings through iLab**, must be made **24 hours in advance** on workdays.
  - **8 hours:** Cartridge loading and SPRINT Run only
  - **24 hours:** Total RNA CodeSet hybridization, cartridge loading, and SPRINT run
  - **24 hours + miRNA Sample Prep:** Contact Jinah about booking extra bench time for the miRNA Sample Prep protocol.
- **Current Run Fee:**
  - \$285 OHSU Internal
  - \$328 External Academic
  - Slated to increase in FY25
- **NanoString Instrument can only be used by trained users.**

# We provide following Support:

- Direct client support for questions and liaison with Nanostring customer service.
- Monitor the instrument during your NanoString run.
- Email Data to the Users.
- Maintain the instrument through the annual service contract provided by NanoString (~\$15,000/year).
- Provide the basic supplies such as calibrated pipets and consumables.



# Support items provided by IGL RNA/DNA Services



- -20 °C storage

# Not provided, must be supplied by user:

- ✓ Samples
- ✓ NanoString CodeSet Panel
- ✓ Sample Prep Kits (if applicable)
- ✓ NanoString SPRINT kit Cartridge
- ✓ Report Library File (RLF) USB drive including your CodeSet
- ✓ -80 °C storage

# Other available Services at IGL to support nCounter SPRINT run

- ☐ Manual and Automated Nucleic Acid Isolations
- ☐ RNA quality assessment with Bioanalyzer
- ☐ DNA quality assessment with Tapestation
- ☐ RNA & DNA UV quantification with NanoDrop
- ☐ RNA & DNA fluorometric quantification with Qubit



# Contact information

## IGL members

- Jinah Kim ([kimjinah@ohsu.edu](mailto:kimjinah@ohsu.edu))
- Britt Daughtry ([daughtry@ohsu.edu](mailto:daughtry@ohsu.edu))
- Robert Searles ([searlesr@ohsu.edu](mailto:searlesr@ohsu.edu))

## NanoString support

- [support.spatial@bruker.com](mailto:support.spatial@bruker.com)
- Westley Heydeck ([Westly.Heydeck@bruker.com](mailto:Westly.Heydeck@bruker.com)): Sr. Field Application Scientist
- Bailey Longoni ([Bailey.Longoni@bruker.com](mailto:Bailey.Longoni@bruker.com)): nCounter Regional Account Manager