
 <p> <b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b>          Technical Supervisor: Andrea E. DeBarber, PhD          Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>          Director/Clinical Consultant: P. Bart Duell, MD          Phone: 503-494-3273       </p>	<p> <b>Sterol Analysis Laboratory</b>          Oregon Health &amp; Science University          3181, SW Sam Jackson Park Road          Portland, OR 97239          Laboratory Phone: 503-494-4593       </p> <p> <b>CAP # 2442607</b> <span style="float: right;"><b>CLIA # 38D06-56829</b></span> </p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

# STEROL ANALYSIS CLINICAL LABORATORY SERVICES GUIDE

For questions concerning laboratory protocols, assay development and validation, quality control, sample handling and assay results, please contact:

Paul Barton Duell, MD  
 Laboratory Director, Sterol Analysis Laboratory  
 Phone: (503) 494-2007

Andrea E. DeBarber, PhD  
 Technical Supervisor, Sterol Analysis Laboratory  
 Phone: (503) 494-4593

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

## TABLE OF CONTENTS

### General Information:


Introduction.....	3
Laboratory Mission.....	3
Laboratory Hours of Operation .....	3
Specimen Collection .....	3
Unacceptable Specimens .....	4
Specimen Labeling, Ordering Tests and Shipping .....	4
Shipping Address.....	5
Assay Timing and Release of Results .....	5
Sample Disposal.....	5
Assay Batching .....	5
Color Coding of Blood Collection Tubes .....	6
Specimen Requirements for Common Tests Requested.....	6

### List of Available Tests:

Plasma/Serum Sterols, Total.....	7
Plasma/Serum 7 $\alpha$ -hydroxy-4-cholesten-3-one, Free .....	8
Plasma/Serum 7 $\alpha$ ,12 $\alpha$ -dihydroxy-4-cholesten-3-one, Free.....	9
Plasma/Serum 5 $\beta$ -cholestane-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,25-tetrol-3-O- $\beta$ -D-glucuronide .....	10
Plasma/Serum cholestane-3 $\beta$ ,5 $\alpha$ ,6 $\beta$ -triol, Free .....	11
Plasma/Serum 7-keto cholesterol, Free.....	12
Urinary 5 $\beta$ -cholestane-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,23S,25-pentol, Total .....	13
Urinary 7-keto-cholesterol, Free.....	14

## GENERAL INFORMATION

### INTRODUCTION

 <p> <b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b>          Technical Supervisor: Andrea E. DeBarber, PhD          Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>          Director/Clinical Consultant: P. Bart Duell, MD          Phone: 503-494-3273       </p>	<p> <b>Sterol Analysis Laboratory</b>          Oregon Health &amp; Science University          3181, SW Sam Jackson Park Road          Portland, OR 97239          Laboratory Phone: 503-494-4593       </p> <p> <b>CAP # 2442607</b> <span style="float: right;"><b>CLIA # 38D06-56829</b></span> </p>
<b>Title: Sterol Analysis Clinical Laboratory Services Guide</b>	

The intention of this manual is to provide guidelines for specimen collection and handling, as well as provide a guide of services currently offered by the Sterol Analysis Laboratory. The Sterol Analysis Laboratory primarily performs quantification of stanols/sterols and intermediates in the bile acid pathway. To request laboratory tests, see section for specimen collection, handling, and storage.

The tests listed in this Sterol Analysis Laboratory Service Guide are only those which are currently performed or have been performed in the past. It is likely that new assays will be developed. Please contact the Laboratory Director or Supervisor for test availability and scheduling.

### STEROL ANALYSIS LABORATORY MISSION

The mission of the Sterol Analysis Laboratory is to:

1. Perform analyses of stanols/sterols and intermediates in the bile acid pathway that are not routinely available in general hospital clinical laboratories or elsewhere in Oregon. We are currently the only laboratory on the West Coast to offer this type of testing.
2. Make laboratory services to measure stanols/sterols and intermediates in the bile acid pathway for diagnostic testing available to all clinicians practicing at OHSU and to those in the surrounding area.
3. Develop new assays for clinical research, which have the potential for wider applicability in diagnosis and patient care.


### STEROL ANALYSIS LABORATORY HOURS OF OPERATION

The hours of operation of the Sterol Analysis Laboratory are 8:00 a.m. to 5:00 p.m. Monday through Friday. **For blood drawn at night or on weekends, it is important to read the sections entitled "Specimen Collection" and "Specimen Labeling, Ordering Tests and Shipping".**

### SPECIMEN COLLECTION

**PLEASE READ CAREFULLY** Blood samples drawn locally after hours should be refrigerated at 4°C and delivered the next working day to Richard Jones Hall, Room 3360 or shipped to the Sterol Analysis Laboratory using overnight delivery. Whole blood specimens should be shipped on ice packs (insulate sample and include one ice-pack, do not freeze). Plasma should be shipped on ice packs or frozen. Urine should be shipped on ice-packs or frozen. Whole blood samples received by the Sterol Analysis Laboratory will be spun and separated. Plasma will be analyzed immediately or frozen prior to analysis.

**NOTE:** All blood, plasma and urine samples should be transported in appropriate biohazard sealed containers that are leak proof.

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

UNACCEPTABLE SPECIMENS (Including, but not limited to, the following):

1. Specimens consisting of citrate, oxalate, or fluoride plasma (blue, gray, or black topped tubes).
2. Specimens without at least two identifiers (i.e. patient name, DOB and/or MRN).
3. Blood/plasma specimens of less than 0.5 mL volume.
4. Specimens that include broken tubes.
5. Specimens not refrigerated or frozen as required.
6. Whole blood specimens drawn more than one day before receipt by the laboratory and/or that are grossly hemolyzed.
7. Specimens that contain any needles or sharps.

This is not intended to imply that all "unacceptable" specimens will be discarded or not analyzed. Requesting health care providers who send unacceptable specimens will be notified no later than the next working day. The phlebotomist, if known, will be notified of the problem.

SPECIMEN LABELING, ORDERING TESTS AND SHIPPING

Samples must be labeled with at least two identifiers, i.e. patient name, DOB and/or MRN.


**NOTE:** ordering using EPIC (OHSU Out-Patient & In-Patient). Currently, we are a reference lab to OHSU Clinical Pathology. Tests performed by the Sterol Analysis Laboratory for OHSU include: plasma  $\beta$ -sitosterol, plasma 7-dehydrocholesterol, plasma 5 $\alpha$ -cholestanol, and urine bile alcohol (5 $\beta$ -cholestane-3 $\alpha$ ,7 $\alpha$ ,12  $\alpha$ ,23S,25-pentol).

The Laboratory Sample Requisition Form (or EPIC order) associated with the sample must provide the following information:

- **Patient Name**
- **Ordering Health Care Provider name**
- **Medical Record Number (MRN)**
- **DOB, Gender**
- **Date of specimen collection**
- **Test requested to be performed**

For samples shipped directly to the Sterol Analysis Laboratory from outside OHSU, please include:

- **A paper Sample Requisition Form, completed and signed by Ordering Health Care Provider**
- **We will need to invoice the referring laboratory or bill the patient directly – if patient is responsible, please provide patient address, phone number and email. To bill to insurance please request OHSU sample requisition form, also found here:**  
<https://www.ohsu.edu/sites/default/files/2023-04/Laboratory%20Requisition%2004-11-23.pdf>.

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

#### SHIPPING ADDRESS:

Attention: Andrea DeBarber (503-494-4593)  
Richard Jones Hall Room 3340/3360, Dock 4, Oregon Health & Science University  
3181 SW Sam Jackson Park Road  
Portland, OR 97239-3098

**NOTE:** Specimens should be shipped by overnight express carrier Monday through Thursday. Saturday, Sunday and holiday deliveries are not possible. Call Laboratory at (503) 494-4593 or email tracking number to [debarber@ohsu.edu](mailto:debarber@ohsu.edu) to notify lab of shipment.

#### ASSAY TIMING AND RELEASE OF RESULTS

Assays performed by the Sterol Analysis Laboratory are analyzed on a batch basis whenever possible. Therefore, turnaround times will vary depending on which assay is requested. Turnaround time for each assay is provided in Laboratory Services Guide. The Sterol Analysis Director or Technical Supervisor/Laboratory Supervisor reviews all results.

If samples are referred to our Laboratory by Clinical Pathology, the final reports are provided to Clinical Pathology and reported in EPIC. For samples submitted directly to the Sterol Analysis Laboratory reports are emailed and/or Faxed to the requesting health care provider and/or referring laboratory.


#### ASSAY BATCHING

The following tests are batched and analyzed by the Sterol Analysis Laboratory for analysis:

Plasma/serum sterols (includes  $\beta$ -sitosterol; 7-dehydrocholesterol,  $5\alpha$ -cholestanol)  
Plasma/serum bile acid pathway intermediates ( $7\alpha,12\alpha$ -dihydroxy-4-cholesten-3-one and  $7\alpha$ -hydroxy-4-cholesten-3-one)  
Plasma/serum bile alcohol ( $5\beta$ -cholestane- $3\alpha,7\alpha,12\alpha,25$ -tetrol-3-O- $\beta$ -D-glucuronide)  
Plasma oxysterols (Cholestane- $3\beta,5\alpha,6\beta$ -triol and 7-Keto Cholesterol)  
Urinary bile alcohol ( $5\beta$ -cholestane- $3\alpha,7\alpha,12\alpha,23S,25$ -pentol)  
Urinary oxysterol (7-Keto Cholesterol)

#### SAMPLE DISPOSAL

The Sterol Analysis Laboratory is faced with space limitations. Therefore, samples are generally not stored for longer than 2 years.

 <b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b> Technical Supervisor: Andrea E. DeBarber, PhD Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a> Director/Clinical Consultant: P. Bart Duell, MD Phone: 503-494-3273	<b>Sterol Analysis Laboratory</b> Oregon Health & Science University 3181, SW Sam Jackson Park Road Portland, OR 97239 Laboratory Phone: 503-494-4593  <b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b>
<b>Title: Sterol Analysis Clinical Laboratory Services Guide</b>	

### COLOR CODING OF BLOOD COLLECTION TUBES


TUBE	ADDITIVE	GENERAL USE
Lavender top	EDTA(K2) Spray Dried/2 or 3 mL	Plasma
Green top	143 IU Sodium Heparin/2 or 3 mL	Plasma
Red top	None	Serum

Prior to use for patient testing the laboratory requires that all consumable material used in testing meets manufacturers or laboratory specifications.

### SPECIMEN REQUIREMENTS FOR COMMON TESTS REQUESTED

(Please read also "Specimen Collection" on page 3)

TEST	MEASURE	TUBE REQUIRED
1. Plasma/serum sterol <sup>1</sup>	$\beta$ -Sitosterol, total	3 mL lavender or green top (plasma) or 3 mL red top (serum)
2. Plasma/serum sterol <sup>1</sup>	7-Dehydrocholesterol, total	3 mL lavender or green top (plasma) or 3 mL red top (serum)
3. Plasma/serum sterol <sup>1</sup>	5 $\alpha$ -Cholesterol, total	3 mL lavender or green top (plasma) or 3 mL red top (serum)
4. Plasma/serum bile acid pathway intermediate <sup>1</sup>	7 $\alpha$ -Hydroxy-4-cholesten-3-one, free	3 mL lavender or green top (plasma) or 3 mL red top (serum)
5. Plasma/serum bile acid pathway intermediate <sup>1</sup>	7 $\alpha$ ,12 $\alpha$ -Dihydroxy-4-cholesten-3-one, free	3 mL lavender or green top (plasma) or 3 mL red top (serum)
6. Plasma/serum bile alcohol <sup>1</sup>	5 $\beta$ -Cholestane-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,25-tetrol-3-O- $\beta$ -D-glucuronide	3 mL lavender or green top (plasma) or 3 mL red top (serum)
7. Plasma/serum oxysterol	Cholestane-3 $\beta$ ,5 $\alpha$ ,6 $\beta$ -triol	3 mL lavender or green top (plasma) or 3 mL red top (serum)
8. Plasma/serum oxysterol	7-Keto-cholesterol	3 mL lavender or green top (plasma) or 3 mL red top (serum)
9. Urinary bile alcohol	5 $\beta$ -Cholestane-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,23S,25-pentol, total	5-10 mL random urine
10. Urinary oxysterol	7-Keto-cholesterol	5-10 mL random urine

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

## LIST OF AVAILABLE TESTS AND METHODOLOGY

### OFFERED THROUGH THE STEROL ANALYSIS CLINICAL LABORATORY:

<b>PLASMA/SERUM STEROLS, TOTAL</b>									
Plasma (preferred)/serum	3 mL lavender or green top (plasma) 3 mL red top (serum) (minimum volume 1 mL)								
<p><b>Method:</b> Hexane extraction of saponified plasma lipids followed by derivatization. Gas chromatography for quantitative determination of elevated 5<math>\alpha</math>-cholestanol (for diagnosis of Cerebrotendinous Xanthomatosis or CTX); 7-dehydrocholesterol (for diagnosis of Smith-Lemli-Opitz Syndrome), <math>\beta</math>-sitosterol (for diagnosis of Sitosterolemia). Quantification is by flame ionization detection or selected-ion monitoring, ion-ratio fragmentometry in the electron impact mode using epicoprostanol as internal standard. Linear up to 4.8 mg/dL without dilution.</p> <p>References: Kelley RI (1995) Clin Chim Acta 236:45, Merkens LS et al (2009) J Pediatr 154:557.</p>									
<p><b>Sample Stability:</b> Plasma/serum is stable up to 17 hours at ambient temperature or refrigerated. Samples can be frozen for long-term storage. Plasma can be shipped overnight on ice-packs or frozen. Whole blood can be shipped overnight on ice-packs (<u>not frozen</u>).</p>									
<p><b>Reference Ranges:</b></p> <table border="0"> <tr> <td>Normal 7-dehydrocholesterol concentration in children age range 0.1-16 years old, n=153</td><td align="right">&lt; 0.10 mg/dL</td></tr> <tr> <td>Normal <math>\beta</math>-sitosterol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369</td><td align="right">0.30 <math>\pm</math> 0.14 mg/dL [0.03-0.94]<sup>1</sup> <b>mean <math>\pm</math> SD [range]</b></td></tr> <tr> <td>Normal 5<math>\alpha</math>-cholestanol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369</td><td align="right">0.27 <math>\pm</math> 0.08 mg/dL [0.10-0.65]<sup>2</sup> <b>mean <math>\pm</math> SD [range]</b></td></tr> <tr> <td>Cholestanol concentration in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36</td><td align="right">2.06 <math>\pm</math> 1.55 mg/dL [0.21-7.10]<sup>3,4</sup> <b>mean [range]</b></td></tr> </table> <p>1. 99<sup>th</sup> percentile concentration of <math>\beta</math>-sitosterol previously reported as 7.5 mg/L (Schaefer et al (2016) Endotext.).  2. 99<sup>th</sup> percentile concentration of 5<math>\alpha</math>-cholestanol previously reported as 5.4 mg/L (Schaefer et al (2016) Endotext.).  3. Note that 5<math>\alpha</math>-cholestanol may be normal in CTX and that medications such as Prednisone, Ezetimibe, and statins may result in normal blood/plasma cholestanol, leading to false negative result for CTX. (Siman-Tov T et al (2006) J Neurol Sci. 243(1-2):83–86, DeBarber AE et al (2024) J Clin Lipidol. 18(3):e465-e476.).  4. Note that cholestanol may also be &gt;1.00 mg/dL for other conditions, for example in cholestasis and liver disease (Koopman BJ et al (1984) Clin Chem Acta 137(3):305-315).</p>		Normal 7-dehydrocholesterol concentration in children age range 0.1-16 years old, n=153	< 0.10 mg/dL	Normal $\beta$ -sitosterol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369	0.30 $\pm$ 0.14 mg/dL [0.03-0.94] <sup>1</sup> <b>mean <math>\pm</math> SD [range]</b>	Normal 5 $\alpha$ -cholestanol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369	0.27 $\pm$ 0.08 mg/dL [0.10-0.65] <sup>2</sup> <b>mean <math>\pm</math> SD [range]</b>	Cholestanol concentration in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36	2.06 $\pm$ 1.55 mg/dL [0.21-7.10] <sup>3,4</sup> <b>mean [range]</b>
Normal 7-dehydrocholesterol concentration in children age range 0.1-16 years old, n=153	< 0.10 mg/dL								
Normal $\beta$ -sitosterol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369	0.30 $\pm$ 0.14 mg/dL [0.03-0.94] <sup>1</sup> <b>mean <math>\pm</math> SD [range]</b>								
Normal 5 $\alpha$ -cholestanol concentration in unaffected individuals age range 0.2-70.3 years old, n= 369	0.27 $\pm$ 0.08 mg/dL [0.10-0.65] <sup>2</sup> <b>mean <math>\pm</math> SD [range]</b>								
Cholestanol concentration in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36	2.06 $\pm$ 1.55 mg/dL [0.21-7.10] <sup>3,4</sup> <b>mean [range]</b>								

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.





Technical Supervisor: Andrea E. DeBarber, PhD  
Email: [debarber@ohsu.edu](mailto:debarber@ohsu.edu)  
Director/Clinical Consultant: P. Bart Duell, MD  
Phone: 503-494-3273

**Sterol Analysis Laboratory**  
Oregon Health & Science University  
3181, SW Sam Jackson Park Road  
Portland, OR 97239  
Laboratory Phone: 503-494-4593

CAP # 2442607

CLIA # 38D06-56829

**Title: Sterol Analysis Clinical Laboratory Services Guide**

**PLASMA/SERUM BILE ACID PATHWAY INTERMEDIATE,  
7 $\alpha$ -HYDROXY-4-CHOLESTEN-3-ONE, FREE**

Plasma (preferred)/serum	3 mL lavender or green top (plasma)
	3 mL red top (serum)
	(minimum volume 0.5 mL)

**Method:** Addition of 7 $\alpha$ -hydroxy-4-cholesten-3-one-d<sub>7</sub> method internal standard in methanol to sample and extracted as previously described. Isotope-dilution quantification (for diagnosis of Cerebrotendinous Xanthomatosis or CTX) is performed using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 7 $\alpha$ C4 (m/z 401.2>177.1 and 401.2>97.0) and internal standard (m/z 408.3>177.1) with a reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6  $\mu$ m) column. Linear up to 1,000 ng/mL without dilution. No significant interference from common over the counter and prescription medications. Carry over response in blanks is < 10% of the LLOQ for the analyte and < 1% for the internal standard

References: DeBarber AE et al (2014) Clin Biochem. 47(9):860-3, Yuan L et al (2016) Bioanalysis. 8:2445-2455, Donato LJ et al (2018) Clin Biochem. 52:106-111.

**Precision:** In-house intra-assay CV < 8% for low QC n=20 samples and < 8% for high QC n=20 samples

**Accuracy:** In-house intra-assay within  $\pm$  10% of nominal value for low and high QC, n=20 samples

**Sample Stability:** Plasma is stable for up to 8 hours at ambient temperature (66 hours refrigerated) or can be frozen for long-term storage. Plasma can be shipped overnight on ice-packs or frozen. Do not keep plasma at room temperature for more than 8 hours.

**Reportable Range:**

7 $\alpha$ -HYDROXY-4-CHOLESTEN-3-ONE	20 ng/mL to 10,000 ng/mL
---------------------------------------	--------------------------

**Reference Ranges:**

Normal concentration in unaffected individuals age range 4.0-70.3 years old, n=14	16.3 $\pm$ 19.5 ng/mL <sup>1</sup> (0.041 $\pm$ 0.049 nmol/mL)
--	---

**mean  $\pm$  SD**

Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36	1,611 ng/mL [27-4,863] (4.022 nmol/mL [0.067-12.139])
--	--

**mean [range]**

1. 7 $\alpha$ -hydroxy-4-cholesten-3-one normal reference range calculated including concentrations <LLOQ. Note that reported normal reference ranges are 22  $\pm$  20 ng/mL (Mignarri A et al (2016) J Inherit Metab Dis 39(1):75-83) and 15  $\pm$  4 ng/mL (Matysik S et al (2011) Chem Phys Lipids 164(6):530-4)

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.





Technical Supervisor: Andrea E. DeBarber, PhD  
Email: [debarber@ohsu.edu](mailto:debarber@ohsu.edu)  
Director/Clinical Consultant: P. Bart Duell, MD  
Phone: 503-494-3273

**Sterol Analysis Laboratory**  
Oregon Health & Science University  
3181, SW Sam Jackson Park Road  
Portland, OR 97239  
Laboratory Phone: 503-494-4593

**CAP # 2442607**

**CLIA # 38D06-56829**

**Title: Sterol Analysis Clinical Laboratory Services Guide**

**PLASMA/SERUM BILE ACID PATHWAY INTERMEDIATE**

**7 $\alpha$ ,12 $\alpha$ -DIHYDROXY-4-CHOLESTEN-3-ONE, FREE**

Plasma (preferred)/serum	3 mL lavender or green top (plasma)
	3 mL red top (serum)
	(minimum volume 0.5 mL)

**Method:** Addition of 7 $\alpha$ ,12 $\alpha$ -dihydroxy-4-cholesten-3-one-d<sub>7</sub> method internal standard in methanol to sample and extracted as previously described. Isotope-dilution quantification (for diagnosis of Cerebrotendinous Xanthomatosis or CTX) is performed using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 7 $\alpha$ 12 $\alpha$ C4 (m/z 417.3>253.1 and 417.3>381.2) and internal standard (m/z 424.3>253.1) with a reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6  $\mu$ m) column. Linear up to 1,000 ng/mL without dilution. No significant interference from common over the counter and prescription medications. Carry over response in blanks is < 6% of the LLOQ for the analyte and < 1% for the internal standard.

References: DeBarber AE et al (2018) J Lipid Res. 59 (11):2214-2222, DeBarber AE et al (2014) J Lipid Res. 55(1):146-54, DeBarber AE et al (2014) Clin Biochem. 47(9):860-3.

**Precision:** In-house intra-assay CV < 7% for low QC n=20 samples and < 7% for high QC n=20 samples

**Accuracy:** In-house intra-assay within  $\pm$  10% of nominal value for low and high QC, n=20 samples

**Sample Stability:** Plasma is stable for up to 8 hours at ambient temperature (66 hours refrigerated) or can be frozen for long-term storage. Plasma can be shipped overnight on ice-packs or frozen. Do not keep plasma at room temperature for more than 8 hours.

**Reportable Range:**

7 $\alpha$ ,12 $\alpha$ -DIHYDROXY-4-CHOLESTEN-3-ONE	10 ng/mL to 10,000 ng/mL
--	--------------------------


**Reference Ranges:**

Normal concentration in unaffected individuals	< 10 ng/mL
age range 4.0-70.3 years old, n=14	(< 0.024 nmol/mL)

Concentrations in CTX-affected untreated individuals	2,340 ng/mL [106-5,632]
age range 0.7-84.5 years old, n=36	(5.617 nmol/mL [0.255-13.519])


**mean [range]**

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.

 <p> <b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b>          Technical Supervisor: Andrea E. DeBarber, PhD          Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>          Director/Clinical Consultant: P. Bart Duell, MD          Phone: 503-494-3273       </p>	<p> <b>Sterol Analysis Laboratory</b>          Oregon Health &amp; Science University          3181, SW Sam Jackson Park Road          Portland, OR 97239          Laboratory Phone: 503-494-4593       </p> <p> <b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b> </p>
<b>Title: Sterol Analysis Clinical Laboratory Services Guide</b>	

<b>PLASMA/SERUM BILE ALCOHOL</b> <b>5<math>\beta</math>-CHOLESTANE-3<math>\alpha</math>,7<math>\alpha</math>,12<math>\alpha</math>,25-TETROL GLUCURONIDE</b>					
Plasma (preferred)/serum	3 m lavender or green top (plasma) 3 mL red top (serum) (minimum volume 0.5 mL)				
<p> <b>Method:</b> Addition of 5<math>\beta</math>-cholestane-3<math>\alpha</math>,7<math>\alpha</math>,12<math>\alpha</math>,25-tetrol-3-O-<math>\beta</math>-D-glucuronide-d<sub>6</sub> method internal standard in methanol to sample and extracted as previously described. Isotope-dilution quantification is performed (for diagnosis of Cerebrotendinous Xanthomatosis or CTX) using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 5<math>\beta</math>-cholestane-3<math>\alpha</math>,7<math>\alpha</math>,12<math>\alpha</math>,25-tetrol-3-O-<math>\beta</math>-D-glucuronide (m/z 611.4&gt;75.0 and 611.4&gt;85.0) and internal standard (m/z 617.4&gt;75.0) with reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6 <math>\mu</math>m) column. Linear up to 5,000 ng/mL without dilution. No significant interference from common over the counter and prescription medications. Carry over response in blanks is &lt; 2% of the LLOQ for the analyte and &lt; 1% for the internal standard.         </p> <p>           Reference: DeBarber AE et al (2018) J Lipid Res. 59 (11):2214-2222.         </p>					
<p> <b>Precision:</b> In-house intra-assay CV &lt; 4% for low QC n=20 samples and &lt; 5% for high QC n=20 samples         </p> <p> <b>Accuracy:</b> In-house intra-assay within <math>\pm</math> 10% of nominal value for low and high QC, n=20 samples         </p> <p> <b>Sample Stability:</b> Plasma is stable for up to 8 hours at ambient temperature (66 hours refrigerated) or can be frozen for long-term storage. Plasma can be shipped overnight on ice-packs or frozen. Do not keep plasma at room temperature for more than 8 hours.         </p>					
<p> <b>Reportable Range:</b> </p> <table> <tr> <td>5<math>\beta</math>-CHOLESTANE-3<math>\alpha</math>,7<math>\alpha</math>,12<math>\alpha</math>,25-TETROL GLUCURONIDE</td><td>100 ng/mL to 50,000 ng/mL</td></tr> </table>		5 $\beta$ -CHOLESTANE-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,25-TETROL GLUCURONIDE	100 ng/mL to 50,000 ng/mL		
5 $\beta$ -CHOLESTANE-3 $\alpha$ ,7 $\alpha$ ,12 $\alpha$ ,25-TETROL GLUCURONIDE	100 ng/mL to 50,000 ng/mL				
<p> <b>Reference Ranges:</b> </p> <table> <tr> <td>Normal concentration in unaffected individuals age range 4.0-70.3 years old, n=14</td><td>           106.4 <math>\pm</math> 4.8 ng/mL            (0.174 <math>\pm</math> 0.008 nmol/mL)  <b>mean <math>\pm</math> SD</b> </td></tr> <tr> <td>Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36</td><td>           7,687 ng/mL [331-27,793]            (12.544 nmol/mL [0.540-45.354])  <b>mean [range]</b> </td></tr> </table>		Normal concentration in unaffected individuals age range 4.0-70.3 years old, n=14	106.4 $\pm$ 4.8 ng/mL (0.174 $\pm$ 0.008 nmol/mL) <b>mean <math>\pm</math> SD</b>	Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36	7,687 ng/mL [331-27,793] (12.544 nmol/mL [0.540-45.354]) <b>mean [range]</b>
Normal concentration in unaffected individuals age range 4.0-70.3 years old, n=14	106.4 $\pm$ 4.8 ng/mL (0.174 $\pm$ 0.008 nmol/mL) <b>mean <math>\pm</math> SD</b>				
Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=36	7,687 ng/mL [331-27,793] (12.544 nmol/mL [0.540-45.354]) <b>mean [range]</b>				

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a> Director/Clinical Consultant: P. Bart Duell, MD Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b> Oregon Health &amp; Science University 3181, SW Sam Jackson Park Road Portland, OR 97239 Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

<p align="center"><b>PLASMA/SERUM OXYSTEROL</b></p> <p align="center"><b>CHOLESTANE-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-TRIOL, FREE</b></p>					
<p>Plasma (preferred)/serum</p>	<p>3 m lavender or green top (plasma) 3 mL red top (serum) (minimum volume 0.5 mL)</p>				
<p><b>Method:</b> Addition of cholestane-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-triol-d<sub>7</sub> method internal standard in methanol to extracted sample followed by derivatization to picolinyl esters [1,2] as previously described. Isotope-dilution quantification is performed (for diagnosis of Nieman Pick Type C Disease or NPC) using positive-mode LC-ESI-MS/MS multiple reaction monitoring of cholestane-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-triol (m/z 631.4&gt;508.4 and 631.4&gt;385.2) and internal standard (m/z 638.4&gt;374.2) with reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6 <math>\mu</math>m) column. Linear up to 700 ng/mL without dilution. No significant interference from co-existing compounds. Carry over response in blanks is &lt; 6% of the LLOQ for the analyte and &lt; 1% for the internal standard.</p> <p>References: 1. Yamamuro D, et al. (2020) J Lipid Res. 61 (9):1287-1299. 2. Chen, L., et al. (2019) Chromatographia. 82, 553–564.</p>					
<p><b>Precision:</b> In-house intra-assay CV &lt; 4.93% for low QC n=20 samples and &lt; 3.22% for high QC n=20 samples</p> <p><b>Accuracy:</b> In-house intra-assay within <math>\pm</math> 14% of nominal value for low and high QC, n=20 samples</p> <p><b>Sample Stability:</b> Plasma/serum should be isolated immediately and is stable for up to 2 hours at ambient temperature or can be frozen for long-term storage. Plasma can be shipped overnight on dry ice, frozen. Do not keep plasma at room temperature for more than 2 hours.</p>					
<p><b>Reportable Range:</b></p> <p align="center">CHOLESTANE-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-TRIOL <span style="float:right">20 ng/mL to 14,000 ng/mL</span></p>					
<p><b>Reference Ranges:</b></p> <table border="0" style="width:100%"> <tr> <td style="width:60%">Normal concentration in unaffected individuals n=20</td><td style="text-align:right"> <math>&lt; 20 \text{ ng/mL}^1</math>  <math>(&lt; 0.048 \text{ nmol/mL})</math> </td></tr> <tr> <td>Reported concentrations in NPC-affected untreated individuals</td><td style="text-align:right"> <math>229 \text{ ng/mL}^2</math>  <math>[24.7 - 489] \text{ and } [21.9 - 963]</math>  <b>mean [range]</b> </td></tr> </table> <p><small>1. cholestane-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-triol reported normal reference ranges are 11.4 - 44.4 ng/mL, mean 29.0 ng/mL (Jiang X et al (2011) J Lipid Res. 52(7):1435-1445) and 3.8 - 39.8 ng/mL (Boenzi S et al (2016) J Lipid Res. 57(3):361-367) 2. cholestane-3<math>\beta</math>,5<math>\alpha</math>,6<math>\beta</math>-triol reported NPC affected reference ranges are 24.7 – 489 ng/mL (Jiang X et al (2011) J Lipid Res. 52(7):1435-1445)) and 21.9 - 963 ng/mL, median 86 ng/mL (P &lt; 0.001) (Boenzi S et al (2016) J Lipid Res. 57(3):361-367)</small></p>		Normal concentration in unaffected individuals n=20	$< 20 \text{ ng/mL}^1$ $(< 0.048 \text{ nmol/mL})$	Reported concentrations in NPC-affected untreated individuals	$229 \text{ ng/mL}^2$ $[24.7 - 489] \text{ and } [21.9 - 963]$ <b>mean [range]</b>
Normal concentration in unaffected individuals n=20	$< 20 \text{ ng/mL}^1$ $(< 0.048 \text{ nmol/mL})$				
Reported concentrations in NPC-affected untreated individuals	$229 \text{ ng/mL}^2$ $[24.7 - 489] \text{ and } [21.9 - 963]$ <b>mean [range]</b>				

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.



Technical Supervisor: Andrea E. DeBarber, PhD  
Email: [debarber@ohsu.edu](mailto:debarber@ohsu.edu)  
Director/Clinical Consultant: P. Bart Duell, MD  
Phone: 503-494-3273

**Sterol Analysis Laboratory**  
Oregon Health & Science University  
3181, SW Sam Jackson Park Road  
Portland, OR 97239  
Laboratory Phone: 503-494-4593

**CAP # 2442607**

**CLIA # 38D06-56829**

**Title: Sterol Analysis Clinical Laboratory Services Guide**

**PLASMA/SERUM OXYSTEROL**

**7-KETO CHOLESTEROL, FREE**

Plasma (preferred)/serum	3 m lavender or green top (plasma)
	3 mL red top (serum)
	(minimum volume 0.5 mL)

**Method:** Addition of 7-keto-cholesterol-d<sub>7</sub> method internal standard in methanol to extracted sample followed by derivatization to picolinyl esters [1,2] as previously described. Isotope-dilution quantification is performed (for diagnosis of Nieman Pick Type C Disease or NPC) using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 7-keto-cholesterol (m/z 506.3>383.3 and 506.3>383.31) and internal standard (m/z 513.4>390.2) with reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6 µm) column. Linear up to 600 ng/mL without dilution. No significant interference from co-existing compounds. Carry over response in blanks is < 9% of the LLOQ for the analyte and < 1% for the internal standard.

**References:**

1. Yamamuro D, et al. (2020) J Lipid Res. 61 (9):1287-1299.
2. Chen, L., et al. (2019) Chromatographia. 82, 553–564.

**Precision:** In-house intra-assay CV < 4.96% for low QC n=20 samples and < 2.37% for high QC n=20 samples

**Accuracy:** In-house intra-assay within ± 13% of nominal value for low and high QC, n=20 samples

**Sample Stability:** Plasma/serum should be isolated immediately and is stable for up to 2 hours at ambient temperature or can be frozen for long-term storage. Plasma can be shipped overnight on dry ice, frozen. Do not keep plasma at room temperature for more than 2 hours.

**Reportable Range:**

7-KETO CHOLESTEROL	0 ng/mL to 11,148 ng/mL
--------------------	-------------------------


**Reference Ranges:**

Normal concentration in unaffected individuals n=20	< 43 ng/mL <sup>1</sup> (< 0.107 nmol/mL)
--	--

Reported concentrations in NPC-affected untreated individuals	80.3 ng/mL <sup>2</sup> [15.1–201] and [16.3–608] <b>mean [range]</b>
---	---


1. 7-Keto-cholesterol reported normal reference ranges are 7.42 – 21.2 ng/mL mean 14.6 ng/mL (Jiang X *et al* (2011) *J Lipid Res.* **52**(7):1435-1445) and 1.1 – 21.9 ng/mL (Boenzi S *et al* (2016) *J Lipid Res.* **57**(3):361-367)  
2. 7-Keto-cholesterol reported NPC affected reference ranges are 15.1–201 ng/mL (Jiang X *et al* (2011) *J Lipid Res.* **52**(7):1435-1445) and 16.3–608 ng/mL, median 55.3 ng/mL (Boenzi S *et al* (2016) *J Lipid Res.* **57**(3):361-367)

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

<p align="center"><b>URINARY BILE ALCOHOL</b></p> <p align="center"><b>5β-CHOLESTANE-3α,7α,12α,23S,25-PENTOL, TOTAL</b></p>					
<p>5-10 mL Random Urine, no preservatives (minimum volume 2 mL)</p>					
<p><b>Method:</b> Addition of 23S-pentol-d6 method internal standard followed by incubation of urine with β-glucuronidase enzyme. Isotope-dilution quantification performed (for diagnosis of Cerebrotendinous Xanthomatosis or CTX) using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 23S-pentol (m/z 453.3&gt;361.4) and internal standard (m/z 459.3&gt;343.5) with a reversed-phase gradient utilizing a 4.6x50mm Biphenyl (2.6 μm) column. Linear up to 5,000 ng/mL without dilution. No significant interference from common over the counter and prescription medications. Carry over response in blanks is &lt; 1% of the LLOQ for the analyte and &lt; 0.5% for the internal standard.</p> <p>References: Shimazu K et al (1986) J Biochem 99: 477, Batta AK et al (1985) J Lipid Res 26: 690.</p>					
<p><b>Precision:</b> In-house intra-assay CV &lt; 4% for low QC n=20 samples and &lt; 3% for high QC n=20 samples</p> <p><b>Accuracy:</b> In-house intra-assay within ± 10% of nominal value for low and high QC, n=20 samples</p> <p><b>Sample Stability:</b> Urine is stable up to 5 days at ambient temperature or refrigerated. Samples can be frozen for long-term storage. Urine can be shipped over-night on ice-packs or frozen.</p>					
<p><b>Reportable Range:</b></p> <table> <tr> <td>5β-CHOLESTANE-3α,7α,12α,23S,25-PENTOL</td><td>200 ng/mL to 250,000 ng/mL</td></tr> </table>		5β-CHOLESTANE-3α,7α,12α,23S,25-PENTOL	200 ng/mL to 250,000 ng/mL		
5β-CHOLESTANE-3α,7α,12α,23S,25-PENTOL	200 ng/mL to 250,000 ng/mL				
<p><b>Reference Ranges:</b></p> <table> <tr> <td>Normal concentration in unaffected individuals, n=20</td><td>&lt; 200 ng/mL</td></tr> <tr> <td>Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=31</td><td>40,135 ng/mL [740-177,506] <b>mean [range]</b></td></tr> </table>		Normal concentration in unaffected individuals, n=20	< 200 ng/mL	Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=31	40,135 ng/mL [740-177,506] <b>mean [range]</b>
Normal concentration in unaffected individuals, n=20	< 200 ng/mL				
Concentrations in CTX-affected untreated individuals age range 0.7-84.5 years old, n=31	40,135 ng/mL [740-177,506] <b>mean [range]</b>				

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.

 <p><b>OREGON HEALTH &amp; SCIENCE UNIVERSITY</b></p> <p>Technical Supervisor: Andrea E. DeBarber, PhD  Email: <a href="mailto:debarber@ohsu.edu">debarber@ohsu.edu</a>  Director/Clinical Consultant: P. Bart Duell, MD  Phone: 503-494-3273</p>	<p><b>Sterol Analysis Laboratory</b>  Oregon Health &amp; Science University  3181, SW Sam Jackson Park Road  Portland, OR 97239  Laboratory Phone: 503-494-4593</p> <p><b>CAP # 2442607</b> <b>CLIA # 38D06-56829</b></p>
<p align="center"><b>Title: Sterol Analysis Clinical Laboratory Services Guide</b></p>	

<p align="center"><b>URINARY OXYSTEROL</b></p> <p align="center"><b>7-KETO CHOLESTEROL, FREE</b></p>			
<p>5-10 mL Random Urine, no preservatives (minimum volume 2 mL)</p>			
<p><b>Method:</b> Addition of 7-keto cholesterol-d<sub>7</sub> method internal standard in methanol to extracted sample followed by derivatization to picolinyl esters [1,2] as previously described. Isotope-dilution quantification is performed using positive-mode LC-ESI-MS/MS multiple reaction monitoring of 7-keto-cholesterol (m/z 506.3&gt;383.3 and 506.3&gt;383.31) and internal standard (m/z 513.4&gt;390.2) with reversed-phase gradient utilizing a 150 x 4.6 mm Biphenyl (2.6 µm) column. Linear up to 60 ng/mL without dilution. No significant interference from co-existing compounds. Carry over response in blanks is &lt; 8% of the LLOQ for the analyte and &lt; 1% for the internal standard.</p> <p>Reference:</p> <ol style="list-style-type: none"> <li>1. Yamamuro D, et al. (2020) J Lipid Res. 61 (9):1287-1299.</li> <li>2. Chen, L., et al. (2019) Chromatographia. 82, 553–564.</li> </ol>			
<p><b>Precision:</b> In-house intra-assay CV &lt; 5.5% for low QC n=20 samples and &lt; 4.86% for high QC n=20 samples</p>			
<p><b>Accuracy:</b> In-house intra-assay within ± 19% of nominal value for low and high QC, n=20 samples</p>			
<p><b>Sample Stability:</b> Urine is stable for up to 2 hours at ambient temperature or can be frozen for long-term storage. Plasma can be shipped overnight on dry ice, frozen. Do not keep urine at room temperature for more than 2 hours.</p>			
<p><b>Reportable Range:</b></p> <table> <tr> <td align="center">7-KETO CHOLESTEROL</td><td align="center">1.6 ng/mL to 12,000 ng/mL</td></tr> </table>		7-KETO CHOLESTEROL	1.6 ng/mL to 12,000 ng/mL
7-KETO CHOLESTEROL	1.6 ng/mL to 12,000 ng/mL		
<p><b>Reference Ranges:</b></p> <table> <tr> <td>Normal concentration in unaffected individuals n=20</td><td align="center">&lt; 1.6 ng/mL</td></tr> </table> <p>Note that no data is available on concentrations in NPC-affected untreated individuals.</p>		Normal concentration in unaffected individuals n=20	< 1.6 ng/mL
Normal concentration in unaffected individuals n=20	< 1.6 ng/mL		

Assay is batched by the Sterol Analysis Laboratory with a resultant turn-around time for reporting of 20 business days.