




1



2



3



4

1. Using the METRx™ System dilators, surgeons gain access to the site for screw placement. This process is completed in two locations. Screws are then placed with the long metal screw extenders.
2. Side view of the CD HORIZON® SEXTANT™ Spinal System's screw extenders and the swinging rod inserter arm. The rod inserter delivers the rod through a small opening in the back and through the heads of the screws implanted in the vertebrae.
3. Final rod and screw implant positions.
4. Three small one-inch scars result from procedures with the CD HORIZON® SEXTANT™ Spinal System.