

# **REQUISITION FORM Telomere Length Measurements**

Today's date:				Store patient sample at room temperature  Do not refrigerate				
PATIENT INFORMATION								
Patient's last name: First: Middle:			Birth Da	ate: / yyyyy	Sex □ M □ F			
Patient ID#:				Sample Collect	Time hh / mm			
		REASO	ON FOR TES	TING				
☐ Bone Marrow Failure	☐ Im	munodeficiency	☐ Lymphoid M	alignancy	☐ Myelo	oid Malignancy		
☐ Pulmonary Fibrosis	☐ Otl	ner Lung Disease	☐ Other, please					
ORDERING INFORMATION								
Physician:				NPI#:				
Hospital:								
Address:								
City:				State: Zip Code:				
The person listed as the Orde	ering Ph	nysician is authorized by law to	o order the test.	Results to be sent b	•			
Authorized Signature (Requ	uired):			Fax:				
				Email:				
			T REQUESTE					
Repeat Diagnostics uses the		SH procedure. To th measurements for total l <b>yn</b>			•	rvice, please contact us.		
l <u> </u>		th measurements for total l <b>yn</b>				T-cells and NK cells		
☐ Medical Consultation		00 for a written evaluation by a						
	inform		clinical history, cur	rrent working diagnos	is, symptoi	ms and lab investigations. If the		
		PATIENT ME	EDICAL INFO	RMATION				
BILLING OPTIONS (We do not invoice healthcare insurance companies)								
Institutional Billing:			Patient Billing Credit card (VISA & MasterCard)					
Hospital:			Name on Credit Card:					
Department:			Address:					
Contact:				City:				
Address:				State: Zip Code:				
City:				Card number:				
State:	State: Zip Code: Exp. Dat		Exp. Date (mmyy):	xp. Date (mmyy): CVC:				
Tel:				Signature of Cardholder:				
Email:				Please charge the above credit card in the amount of \$				



## TELOMERE LENGTH MEASUREMENTS SPECIMEN COLLECTION AND SHIPPING PROCEDURE

#### **BEFORE COLLECTION OF BLOOD**

ample should <b>only</b> be collected and shipped on Monday, Tuesday or Wednesday.	
equisition Form check list	
☐ Patient name is filled in and matches blood tube ID (first identifier)	
☐ Second patient identifier (date of birth, unique ID number) is filled in and matches blood tube	
☐ Ordering information is complete and signed by the requesting physician	
Result send out information is completed	
☐ Assay type (2 or 6-panel) and optional consultation are selected accordingly	
☐ Payment information is completed	

#### SPECIMEN COLLECTION

Label the specimen tube with:
□ Patient Name and ID #
□ Age
□ Sex
Date and time of collection
Collect blood in EDTA anti-coagulant tube.
5-10ml of blood is required for successful testing.
Store patient sample at <b>room temperature</b> until pick-up by courier.
All blood shipments to Repeat Diagnostics must arrive within 2 days and in good condition.



#### SPECIMEN PACKING AND SHIPPING

### SHIPPING MATERIAL

UN3373 shipping box measuring approximately 9" X 4" X 4", labeled "E	Biological Substance Category B)
Specimen bag or sealable plastic bag.	
Absorbent material such as paper towel.	*
Packing tape.	COURT SUBSTANCE, CATEGORY B
Address label.	TC-125-18
FedEx Clinical Pak (provided free of charge from FedEx)	UN3373/
International Air Waybill.	100
Commercial Invoice.	
For more information on how to ship clinical samples visit FedEx at http://images.fedex.com/downloads/shared/packagingtips/pointers.pd	lf .

## **SHIPPING**

- 1. Place blood collection tube(s) in sealable plastic bag.
- Place bag in shipping container. ICE PACKS ARE NOT REQUIRED 2.
- 3. Place enough absorbent material in shipping container so that blood tubes do not roll around.
- Seal shipping container with packing tape. 4.
- 5. Attach address label to top of shipping container.
- 6. Place shipping container and requisition form inside FedEx Clinical Pak.
- 7. Fill out the international Air Waybill form.
- 8. Fill out commercial invoice form. Minimal dollar value must be \$4.00 to ensure rapid customs processing.
- Include 5 copies of the Commercial Invoice with the waybill. 9.
- Ship on day of collection by FedEx International Priority to:

Repeat Diagnostics Inc. Suite 309 - 267 West Esplanade North Vancouver, BC V7M 1A5

11. Inform Repeat Diagnostics by email at test@repeatdx.com of date shipped and tracking number.

			COMMERO	CIALIN	v o	ICE			
Date of Export	tation:			Export References :		Clinical Diag	nostic Test		
Shipper/Exporter (complete name and address)			Consignee:  Repeat Diagnostics Suite 309 267 West Esplanade North Vancouver, BC V7M 1A5 Canada  T. 604-985-2609 F. 778-340-1144						
Country of Export United States			Importer same as consignee: Repeat Diagnostics Customs Broker is:						
Country of Ori	gin of Good	s United	States	FedEx EXPRESSCLEAR  Vancouver BC Canada					
Country of Ultin	mate Destina	ation Can	ada						
			International Air Waybill No.						
Marks/Nos.	No. of Pkgs.	Type of Packaging	Full Description of Good	ds	Qty.	Unit of Measure	Weight	Unit Value	Total Value
	1	Box	Fresh Cells Human White Blood C	ells	1		0.5 kilo	4.00	4.00
			For Diagnostics Testing						0.00
			Non-infectious/Non-hazardous/Non	-toxic/Non-volatile					0.00
			No Commercial Value						0.00
									0.00
									0.00
TOTAL	1						0.5 Kgs		\$4.00
Diversion contra	ary to United	States law is pro	nate Destination shown.  hibited.  invoice to be true and correct.			Check On	e	F.O.B.  C&F  C.I.F.	
	Signature	of Shipper							
	(Type nan	ne and title)				Date			